

Product datasheet for **TR300319**

ZFYVE19 Human shRNA Plasmid Kit (Locus ID 84936)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ZFYVE19 Human shRNA Plasmid Kit (Locus ID 84936)
Locus ID:	84936
Synonyms:	ANCHR; MPFYVE
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ZFYVE19 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 84936). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC015738 , BC021092 , NM_001077268 , NM_001258420 , NM_001258421 , NM_032850 , NM_001077268.1 , NM_001258421.1 , NM_001258420.1 , NM_032850.1 , NM_032850.2 , NM_032850.3 , NM_032850.4 , BC021092.1 , NM_001258420.2 , NM_001077268.2 , NM_032850.5 , NM_001258421.2
UniProt ID:	Q96K21
Summary:	Key regulator of abscission step in cytokinesis: part of the cytokinesis checkpoint, a process required to delay abscission to prevent both premature resolution of intercellular chromosome bridges and accumulation of DNA damage. Together with CHMP4C, required to retain abscission-competent VPS4 (VPS4A and/or VPS4B) at the midbody ring until abscission checkpoint signaling is terminated at late cytokinesis. Deactivation of AURKB results in dephosphorylation of CHMP4C followed by its dissociation from ZFYVE19/ANCHR and VPS4 and subsequent abscission.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).