

## Product datasheet for **TR30035V**

### Lentiviral Control Particles


#### Product data:

Product Type: Lentiviral Control Particles

Description: Lenti scrambled control particles of pmCherry-CB-shlenti vector, >1x10<sup>7</sup> TU/ml, 0.5 ml.

Fluorescent Label: None

#### Product images:



mCherry image (10x). High mCherry expression levels and high transduction efficiency were observed. (Criteria: >50% RFP+ in HEK293T cells; Result: >80% RFP+, Pass)



[View online »](#)