

## Product datasheet for **TR30035V5**

### Lentiviral Control Particles


#### Product data:

Product Type: Lentiviral Control Particles

Description: Lenti scrambled control particles of pmCherry-CB-shlenti vector,  $>1 \times 10^7 \mu\text{TU/ml}$ , 0.5 ml x 5.

Fluorescent Label: None

#### Product images:



mCherry image (10x). High mCherry expression levels and high transduction efficiency were observed. (Criteria:  $>50\%$  RFP+ in HEK293T cells; Result:  $>80\%$  RFP+, Pass)



[View online »](#)