


Product datasheet for TR30036

Lentiviral Packaging Kits

Product data:

Product Type:	Lentiviral Packaging Kits
Description:	IDLV Lenti packaging kit, integration deficient- packaging plasmids and transfection reagent
Components:	1 vial of packaging plasmid TR30036A (60 ug) – lyophilized 1 vial of TurboFectin transfection reagent (400 ul) – liquid
Storage:	Upon receipt, store the packaging plasmids and transfection reagent at 4°C.
Stability:	Stable for 1 year from receipt of products under proper storage and handling condition.

Product images:



GFP lentivirus transduced HEK293 cells, normal lentivirus (WT) or integration-deficient (IDLV). GFP positive cells were measured with FACS sorting.



[View online »](#)

Effective KO of RFP W/O Footprint of Cas9

