

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


# Product datasheet for TR30036P5

### **Lentiviral Packaging Kits**

#### **Product data:**

Product Type: Description:	Lentiviral Packaging Kits IDLV Lenti packaging kit, integration deficient- packaging plasmids and transfection reagent, 5
Components:	packs Five kits each containing: 1 vial of packaging plasmid TR30036A (60 ug) – lyophilized 1 vial of Turbofectin (400 ul) – liquid
Storage:	Store in 4°C refrigerator upon arrival.
Stability:	Stable for 1 year from receipt of products under proper storage and handling condition.


## **Product images:**



GFP lentivirus transduced HEK293 cells, normal lentivirus (WT) or integration-deficient (IDLV). GFP positive cells were measured with FACS sorting.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US