


Product datasheet for **TR30037**

Lentiviral Packaging Kits

Product data:

Product Type:	Lentiviral Packaging Kits
Description:	Lenti-vpak packaging kit - packaging plasmids and transfection reagent
Components:	1 vial of packaging plasmids (60 ug) – lyophilized 1 vial of TurboFectin transfection reagent (400 ul) – liquid
Storage:	Upon receipt, store the packaging plasmids and transfection reagent at 4°C.
Stability:	Stable for 1 year from receipt of products under proper storage and handling condition.

Product images:



[View online »](#)