

Product datasheet for **TR30040P5**

Lentiviral Packaging Kits

Product data:

Product Type:	Lentiviral Packaging Kits
Description:	Lenti-vpak packaging kit - packaging plasmids and transfection reagent, 5 packs
Components:	Five kits each containing: 1 vial of packaging plasmid TR30037A (60 ug) – lyophilized 1 vial of Megatran transfection reagent (400 ul) – liquid
Storage:	Store in 4°C refrigerator upon arrival.
Stability:	Stable for 1 year from receipt of products under proper storage and handling condition.



[View online »](#)