

## Product datasheet for **TR300509**

### WDR4 Human shRNA Plasmid Kit (Locus ID 10785)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	WDR4 Human shRNA Plasmid Kit (Locus ID 10785)
Locus ID:	10785
Synonyms:	GAMOS6; hWH; MIGSB; TRM82; TRMT82; Wuho
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	WDR4 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 10785). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001260474</a> , <a href="#">NM_001260475</a> , <a href="#">NM_001260476</a> , <a href="#">NM_001260477</a> , <a href="#">NM_018669</a> , <a href="#">NM_033661</a> , <a href="#">NM_033662</a> , <a href="#">NR_048535</a> , <a href="#">NM_018669.1</a> , <a href="#">NM_018669.2</a> , <a href="#">NM_018669.3</a> , <a href="#">NM_018669.4</a> , <a href="#">NM_018669.5</a> , <a href="#">NM_033661.1</a> , <a href="#">NM_033661.2</a> , <a href="#">NM_033661.3</a> , <a href="#">NM_033661.4</a> , <a href="#">NM_001260475.1</a> , <a href="#">NM_001260476.1</a> , <a href="#">NM_001260477.1</a> , <a href="#">NM_001260474.1</a> , <a href="#">BC006341</a> , <a href="#">BC006341.2</a> , <a href="#">BC001074</a> , <a href="#">BC018866</a> , <a href="#">NM_018669.6</a>
UniProt ID:	<a href="#">P57081</a>
Summary:	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartate (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This gene is excluded as a candidate for a form of nonsyndromic deafness (DFNB10), but is still a candidate for other disorders mapped to 21q22.3 as well as for the development of Down syndrome phenotypes. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, May 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online >](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).