

Product datasheet for **TR300758**

TTC26 Human shRNA Plasmid Kit (Locus ID 79989)

Product data:

Product Type:	shRNA Plasmids
Product Name:	TTC26 Human shRNA Plasmid Kit (Locus ID 79989)
Locus ID:	79989
Synonyms:	dyf-13; DYF13; IFT56
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TTC26 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 79989). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001144920 , NM_001144923 , NM_001287512 , NM_001287513 , NM_001318333 , NM_001321740 , NM_001321741 , NM_001321742 , NM_024926 , NM_024926.1 , NM_024926.2 , NM_024926.3 , NM_001144920.1 , NM_001144920.2 , NM_001144923.1 , NM_001144923.2 , NM_001287512.1 , NM_001287513.1 , BC130339 , BC013912 , BC034466 , BC126331 , BC144151 , NM_001287513.2 , NM_001144923.3 , NM_001287512.2 , NM_024926.4
UniProt ID:	A0AVF1
Summary:	Component of the intraflagellar transport (IFT) complex B required for transport of proteins in the motile cilium. Required for transport of specific ciliary cargo proteins related to motility, while it is neither required for IFT complex B assembly or motion nor for cilium assembly. Required for efficient coupling between the accumulation of GLI2 and GLI3 at the ciliary tips and their dissociation from the negative regulator SUFU. Plays a key role in maintaining the integrity of the IFT complex B and the proper ciliary localization of the IFT complex B components. Not required for IFT complex A ciliary localization or function. Essential for maintaining proper microtubule organization within the ciliary axoneme.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).