

Product datasheet for **TR301186**

TBXAS1 Human shRNA Plasmid Kit (Locus ID 6916)

Product data:

Product Type:	shRNA Plasmids
Product Name:	TBXAS1 Human shRNA Plasmid Kit (Locus ID 6916)
Locus ID:	6916
Synonyms:	BDPLT14; CYP5; CYP5A1; GHOSAL; THAS; TS; TXAS; TXS
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TBXAS1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6916). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001061 , NM_001130966 , NM_001166253 , NM_001166254 , NM_001314028 , NM_030984 , NR_029394 , NM_001061.1 , NM_001061.2 , NM_001061.3 , NM_001061.4 , NM_030984.1 , NM_030984.2 , NM_030984.3 , NM_001130966.1 , NM_001130966.2 , NM_001166253.1 , NM_001166254.1 , BC041157 , BC041157.1 , BC014117 , BM673654 , NM_001366538 , NR_159363 , NR_159365 , NM_001366537 , NR_159364 , NM_001166254.3
UniProt ID:	P24557
Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. However, this protein is considered a member of the cytochrome P450 superfamily on the basis of sequence similarity rather than functional similarity. This endoplasmic reticulum membrane protein catalyzes the conversion of prostglandin H2 to thromboxane A2, a potent vasoconstrictor and inducer of platelet aggregation. The enzyme plays a role in several pathophysiological processes including hemostasis, cardiovascular disease, and stroke. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).