

Product datasheet for **TR301367**

ST3GAL3 Human shRNA Plasmid Kit (Locus ID 6487)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ST3GAL3 Human shRNA Plasmid Kit (Locus ID 6487)
Locus ID:	6487
Synonyms:	DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ST3GAL3 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6487). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001270459 , NM_001270460 , NM_001270461 , NM_001270462 , NM_001270463 , NM_001270464 , NM_001270465 , NM_001270466 , NM_006279 , NM_174963 , NM_174964 , NM_174965 , NM_174966 , NM_174967 , NM_174968 , NM_174969 , NM_174970 , NM_174971 , NM_174972 , NR_073016 , NR_073017 , NR_073018 , NR_073019 , NR_073020 , NR_073021 , NR_073023 , NM_001350619 , NM_001350620 , NM_001350621 , NR_146867 , NM_006279.1 , NM_006279.2 , NM_006279.3 , NM_174966.1 , NM_174966.2 , NM_174963.1 , NM_174963.2 , NM_174963.3 , NM_174969.1 , NM_174969.2 , NM_174964.1 , NM_174964.2 , NM_174965.1 , NM_174965.2 , NM_174967.1 , NM_174967.2 , NM_174968.1 , NM_174968.2 , NM_174968.3 , NM_174970.1 , NM_174970.2 , NM_174971.1 , NM_174971.2 , NM_174971.3 , NM_001270466.1 , NM_001270465.1 , NM_001270464.1 , NM_001270463.1 , NM_001270462.1 , NM_001270461.1 , NM_001270460.1 , NM_001270459.1 , BC050380 , BC050380.1 , NM_001363573 , NM_174967.4 , NM_001270461.3 , NM_174966.4 , NM_174970.4 , NM_001270460.2 , NM_001270462.3 , NM_001270464.3 , NM_174969.4 , NM_174965.3 , NM_006279.5 , NM_174964.4 , NM_001270459.2 , NM_001270466.2 , NM_174963.5 , NM_001270463.3 , NM_001270465.2 , NM_174968.5 , NM_174971.5
UniProt ID:	Q11203



[View online »](#)

Summary:	<p>The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Mutations in this gene have been associated with a form of autosomal recessive nonsyndromic cognitive disability as well as infantile epileptic encephalopathy. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>