

## Product datasheet for **TR301442**

### SPAG6 Human shRNA Plasmid Kit (Locus ID 9576)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SPAG6 Human shRNA Plasmid Kit (Locus ID 9576)
Locus ID:	9576
Synonyms:	CFAP194; CT141; FAP194; pf16; Repro-SA-1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	SPAG6 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9576). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001253854</a> , <a href="#">NM_001253855</a> , <a href="#">NM_012443</a> , <a href="#">NM_172242</a> , <a href="#">NM_172242.1</a> , <a href="#">NM_172242.2</a> , <a href="#">NM_012443.1</a> , <a href="#">NM_012443.2</a> , <a href="#">NM_012443.3</a> , <a href="#">NM_001253854.1</a> , <a href="#">NM_001253855.1</a> , <a href="#">BC030585</a> , <a href="#">BC030585.2</a> , <a href="#">NM_001253854.2</a> , <a href="#">NM_012443.4</a> , <a href="#">NM_172242.3</a>
UniProt ID:	<a href="#">O75602</a>
Summary:	The correlation of anti-sperm antibodies with cases of unexplained infertility implicates a role for these antibodies in blocking fertilization. Improved diagnosis and treatment of immunologic infertility, as well as identification of proteins for targeted contraception, are dependent on the identification and characterization of relevant sperm antigens. The protein expressed by this gene is recognized by anti-sperm antibodies from an infertile man. This protein localizes to the tail of permeabilized human sperm and contains eight contiguous armadillo repeats, a motif known to mediate protein-protein interactions. Studies in mice suggest that this protein is involved in sperm flagellar motility and maintenance of the structural integrity of mature sperm. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).