

Product datasheet for **TR301964**

RNF14 Human shRNA Plasmid Kit (Locus ID 9604)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RNF14 Human shRNA Plasmid Kit (Locus ID 9604)
Locus ID:	9604
Synonyms:	ARA54; HFB30; HRIHFB2038; TRIAD2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RNF14 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9604). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001201365 , NM_004290 , NM_183398 , NM_183399 , NM_183400 , NM_183401 , NM_004290.1 , NM_004290.2 , NM_004290.4 , NM_183400.1 , NM_183400.2 , NM_183401.1 , NM_183401.2 , NM_183398.1 , NM_183398.2 , NM_183399.1 , NM_183399.2 , NM_001201365.1 , BC126185 , BC144061 , NM_183400.3 , NM_183401.3 , NM_004290.5
UniProt ID:	Q9UBS8
Summary:	The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. This protein interacts with androgen receptor (AR) and may function as a coactivator that induces AR target gene expression in prostate. A dominant negative mutant of this gene has been demonstrated to inhibit the AR-mediated growth of prostate cancer. This protein also interacts with class III ubiquitin-conjugating enzymes (E2s) and may act as a ubiquitin-ligase (E3) in the ubiquitination of certain nuclear proteins. Six alternatively spliced transcript variants encoding two distinct isoforms have been reported. [provided by RefSeq, Jan 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).