

Product datasheet for **TR302111**

RAPH1 Human shRNA Plasmid Kit (Locus ID 65059)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RAPH1 Human shRNA Plasmid Kit (Locus ID 65059)
Locus ID:	65059
Synonyms:	LPD, RMO1, PREL2, ALS2CR9, ALS2CR18, KIAA1681, RaIGDS/AF-6
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RAPH1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 65059). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_025252 , NM_203365 , NM_213589 , NM_001329728 , NM_203365.1 , NM_203365.2 , NM_203365.3 , NM_213589.1 , NM_213589.2 , NM_025252.3 , BC156182 , BC156922 , NM_213589.3 , NM_203365.4
UniProt ID:	Q70E73
Summary:	This gene encodes a protein that belongs to the Mig10/Rap1-interacting adaptor molecule/Lamellipodin family of adapter proteins, which function in cell migration. Members of this family contain pleckstrin-homology domains, Ras-association domains, and proline-rich C-termini. The protein encoded by this gene regulates actin dynamics through interaction with Ena/Vasodilator proteins as well as direct binding to filamentous actin to regulate actin network assembly. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).