

## Product datasheet for **TR302358**

### PPAN Human shRNA Plasmid Kit (Locus ID 56342)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PPAN Human shRNA Plasmid Kit (Locus ID 56342)
Locus ID:	56342
Synonyms:	BXDC3; SSF; SSF-1; SSF1; SSF2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	PPAN - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 56342). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_020230</a> , <a href="#">NM_001346139</a> , <a href="#">NM_001346140</a> , <a href="#">NM_001346141</a> , <a href="#">NM_020230.1</a> , <a href="#">NM_020230.2</a> , <a href="#">NM_020230.3</a> , <a href="#">NM_020230.4</a> , <a href="#">NM_020230.5</a> , <a href="#">NM_020230.6</a> , <a href="#">BC033202</a> , <a href="#">BC033202.1</a> , <a href="#">BC000535</a> , <a href="#">BC009833</a> , <a href="#">NM_020230.7</a>
UniProt ID:	<a href="#">Q9NQ55</a>
Summary:	The protein encoded by this gene is an evolutionarily conserved protein similar to yeast SSF1 as well as to the gene product of the Drosophila gene peter pan (ppan). SSF1 is known to be involved in the second step of mRNA splicing. Both SSF1 and ppan are essential for cell growth and proliferation. Exogenous expression of this gene was reported to reduce the anchorage-independent growth of some tumor cells. Read-through transcription of this gene with P2RY11/P2Y(11), an adjacent downstream gene that encodes an ATP receptor, has been found. These read-through transcripts are ubiquitously present and up-regulated during granulocyte differentiation. [provided by RefSeq, Nov 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).