

Product datasheet for **TR302415**

PLP1 Human shRNA Plasmid Kit (Locus ID 5354)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PLP1 Human shRNA Plasmid Kit (Locus ID 5354)
Locus ID:	5354
Synonyms:	GPM6C; HLD1; MMPL; PLP; PLP/DM20; PMD; SPG2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	PLP1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 5354). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_000533 , NM_001128834 , NM_001305004 , NM_199478 , NM_000533.1 , NM_000533.2 , NM_000533.3 , NM_000533.4 , NM_199478.1 , NM_199478.2 , NM_001128834.1 , NM_001128834.2 , BC002665 , BC095452 , NM_199478.3 , NM_000533.5
UniProt ID:	P60201
Summary:	This gene encodes a transmembrane proteolipid protein that is the predominant component of myelin. The encoded protein may play a role in the compaction, stabilization, and maintenance of myelin sheaths, as well as in oligodendrocyte development and axonal survival. Mutations in this gene cause Pelizaeus-Merzbacher disease and spastic paraplegia type 2. Alternatively splicing results in multiple transcript variants, including the DM20 splice variant. [provided by RefSeq, Feb 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).