

## Product datasheet for **TR302418**

### PLEKHO1 Human shRNA Plasmid Kit (Locus ID 51177)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	PLEKHO1 Human shRNA Plasmid Kit (Locus ID 51177)
Locus ID:	51177
Synonyms:	CKIP-1; CKIP1; JBP; OC120
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	PLEKHO1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 51177). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001304722</a> , <a href="#">NM_001304723</a> , <a href="#">NM_001304724</a> , <a href="#">NM_016274</a> , <a href="#">NM_016274.1</a> , <a href="#">NM_016274.2</a> , <a href="#">NM_016274.3</a> , <a href="#">NM_016274.4</a> , <a href="#">NM_016274.5</a> , <a href="#">BC010149</a> , <a href="#">BC010149.2</a> , <a href="#">BC023533</a> , <a href="#">NM_016274.6</a>
UniProt ID:	<a href="#">Q53GL0</a>
Summary:	Plays a role in the regulation of the actin cytoskeleton through its interactions with actin capping protein (CP). May function to target CK2 to the plasma membrane thereby serving as an adapter to facilitate the phosphorylation of CP by protein kinase 2 (CK2). Appears to target ATM to the plasma membrane. Appears to also inhibit tumor cell growth by inhibiting AKT-mediated cell-survival. Also implicated in PI3K-regulated muscle differentiation, the regulation of AP-1 activity (plasma membrane bound AP-1 regulator that translocates to the nucleus) and the promotion of apoptosis induced by tumor necrosis factor TNF. When bound to PKB, it inhibits it probably by decreasing PKB level of phosphorylation.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).