

Product datasheet for **TR302505**

PHF7 Human shRNA Plasmid Kit (Locus ID 51533)

Product data:

Product Type:	shRNA Plasmids
Product Name:	PHF7 Human shRNA Plasmid Kit (Locus ID 51533)
Locus ID:	51533
Synonyms:	HSPC045; HSPC226; NYD-SP6
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	PHF7 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 51533). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001278221 , NM_001321126 , NM_001321127 , NM_016483 , NM_173341 , NM_016483.1 , NM_016483.2 , NM_016483.3 , NM_016483.4 , NM_016483.5 , NM_016483.6 , NM_173341.1 , NM_001278221.1 , NM_001278221.2 , BC022002 , BC022002.1 , NM_001278221.3 , NM_016483.7
UniProt ID:	Q9BWX1
Summary:	Spermatogenesis is a complex process regulated by extracellular and intracellular factors as well as cellular interactions among interstitial cells of the testis, Sertoli cells, and germ cells. This gene is expressed in the testis in Sertoli cells but not germ cells. The protein encoded by this gene contains plant homeodomain (PHD) finger domains, also known as leukemia associated protein (LAP) domains, believed to be involved in transcriptional regulation. The protein, which localizes to the nucleus of transfected cells, has been implicated in the transcriptional regulation of spermatogenesis. Alternate splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).