

## Product datasheet for **TR303682**

### KIRREL2 Human shRNA Plasmid Kit (Locus ID 84063)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	KIRREL2 Human shRNA Plasmid Kit (Locus ID 84063)
Locus ID:	84063
Synonyms:	FILTRIN; NEPH3; NLG1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	KIRREL2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 84063). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC064925</a> , <a href="#">NM_001329530</a> , <a href="#">NM_032123</a> , <a href="#">NM_199179</a> , <a href="#">NM_199180</a> , <a href="#">NM_032123.1</a> , <a href="#">NM_032123.2</a> , <a href="#">NM_032123.3</a> , <a href="#">NM_032123.5</a> , <a href="#">NM_032123.6</a> , <a href="#">NM_199179.1</a> , <a href="#">NM_199179.2</a> , <a href="#">NM_199179.3</a> , <a href="#">NM_199180.1</a> , <a href="#">NM_199180.2</a> , <a href="#">BC007312</a> , <a href="#">BM693057</a> , <a href="#">NM_001363667</a> , <a href="#">NM_199180.4</a> , <a href="#">NM_032123.7</a> , <a href="#">NM_199179.4</a>
UniProt ID:	<a href="#">Q6UWL6</a>
Summary:	This gene encodes a type I transmembrane protein and member of the immunoglobulin superfamily of cell adhesion molecules. The encoded protein localizes to adherens junctions in pancreatic beta cells and regulates insulin secretion. Autoantibodies against the encoded protein have been detected in serum from patients with type 1 diabetes. This gene may also play a role in glomerular development and decreased expression of this gene has been observed in human glomerular diseases. This gene and the related opposite-strand gene nephrin (GeneID: 527362) are regulated by a bidirectional promoter. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).