

## Product datasheet for **TR303835**

### **Kv beta 1 (KCNAB1) Human shRNA Plasmid Kit (Locus ID 7881)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	Kv beta 1 (KCNAB1) Human shRNA Plasmid Kit (Locus ID 7881)
Locus ID:	7881
Synonyms:	AKR6A3; hKvb3; hKvBeta3; KCNA1B; KV-BETA-1; Kvb1.3
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	KCNAB1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 7881). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001308217</a> , <a href="#">NM_001308222</a> , <a href="#">NM_003471</a> , <a href="#">NM_172159</a> , <a href="#">NM_172160</a> , <a href="#">NM_172159.1</a> , <a href="#">NM_172159.2</a> , <a href="#">NM_172159.3</a> , <a href="#">NM_003471.1</a> , <a href="#">NM_003471.2</a> , <a href="#">NM_003471.3</a> , <a href="#">NM_172160.1</a> , <a href="#">NM_172160.2</a> , <a href="#">BC043166</a> , <a href="#">BC043166.2</a> , <a href="#">NM_172160.3</a>
UniProt ID:	<a href="#">Q14722</a>
Summary:	Potassium channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in Drosophila, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shaker-related subfamily. This member includes distinct isoforms which are encoded by alternatively spliced transcript variants of this gene. Some of these isoforms are beta subunits, which form heteromultimeric complexes with alpha subunits and modulate the activity of the pore-forming alpha subunits. [provided by RefSeq, Apr 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).