

Product datasheet for **TR304008**

HVCN1 Human shRNA Plasmid Kit (Locus ID 84329)

Product data:

Product Type:	shRNA Plasmids
Product Name:	HVCN1 Human shRNA Plasmid Kit (Locus ID 84329)
Locus ID:	84329
Synonyms:	HV1; VSOP
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	HVCN1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 84329). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC009731 , NM_001040107 , NM_001256413 , NM_032369 , NM_001040107.1 , NM_032369.1 , NM_032369.2 , NM_032369.3 , NM_001256413.1 , BC009731.1 , BC007277 , BC032672 , NM_001040107.2 , NM_032369.4 , NM_001256413.2
UniProt ID:	Q96D96
Summary:	This gene encodes a voltage-gated protein channel protein expressed more highly in certain cells of the immune system. Phagocytic cells produce superoxide anions which require this channel protein, and in B cells this same process facilitates antibody production. This same channel protein, however, can also regulate functions in other cells including spermatozoa. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).