

## Product datasheet for **TR304126**

### HERC4 Human shRNA Plasmid Kit (Locus ID 26091)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	HERC4 Human shRNA Plasmid Kit (Locus ID 26091)
Locus ID:	26091
Synonyms:	KIAA1593, DKFZP564G092
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	HERC4 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 26091). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC015043</a> , <a href="#">BC039600</a> , <a href="#">NM_001017972</a> , <a href="#">NM_001278185</a> , <a href="#">NM_001278186</a> , <a href="#">NM_001278187</a> , <a href="#">NM_015601</a> , <a href="#">NM_022079</a> , <a href="#">NM_015601.1</a> , <a href="#">NM_015601.2</a> , <a href="#">NM_015601.3</a> , <a href="#">NM_022079.1</a> , <a href="#">NM_022079.2</a> , <a href="#">NM_001278187.1</a> , <a href="#">NM_001278186.1</a> , <a href="#">NM_001278185.1</a> , <a href="#">BC039600.1</a> , <a href="#">NM_001017972.1</a> , <a href="#">BC018114</a> , <a href="#">BC039465</a> , <a href="#">NM_001278187.2</a> , <a href="#">NM_015601.4</a> , <a href="#">NM_001278185.2</a> , <a href="#">NM_001278186.2</a> , <a href="#">NM_022079.3</a>
UniProt ID:	<a href="#">Q5GLZ8</a>
Summary:	HERC4 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).