

Product datasheet for TR304489

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335)

Product data:

Product Type: shRNA Plasmids

Product Name: Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335)

Locus ID: 2335

Synonyms: CIG; ED-B; FINC; FN; FNZ; GFND; GFND2; LETS; MSF; SMDCF

Vector: pRS (TR20003)

E. coli Selection: Ampicillin

Mammalian Cell Puromycin

Selection:

Format: Retroviral plasmids

Components: FN1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID =

2335). 5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.

RefSeq: NM 001306129, NM 001306130, NM 001306131, NM 001306132, NM 002026, NM 054034,

NM 212474, NM 212475, NM 212476, NM 212478, NM 212482, NM 054034.1, NM 054034.2, NM 212482.1, NM 212482.2, NM 212474.1, NM 212474.2, NM 212476.1, NM 212476.2, NM 212478.1, NM 212478.2, NM 002026.1, NM 002026.2, NM 002026.3, NM 212475.1, BC117176, BC005858, BC016875, BC078656, BC100030, BC143754, BC143763, NM 001365517,

NM 001365520, NM 001365518, NM 001365519, NM 001365521, NM 001365522, NM 001365523, NM 001365524, NM 212482.3, NM 054034.3, NM 212474.3

UniProt ID: P02751

Summary: This gene encodes fibronectin, a glycoprotein present in a soluble dimeric form in plasma,

and in a dimeric or multimeric form at the cell surface and in extracellular matrix. The encoded preproprotein is proteolytically processed to generate the mature protein. Fibronectin is involved in cell adhesion and migration processes including embryogenesis, wound healing, blood coagulation, host defense, and metastasis. The gene has three regions

subject to alternative splicing, with the potential to produce 20 different transcript variants, at least one of which encodes an isoform that undergoes proteolytic processing. The full-length

nature of some variants has not been determined. [provided by RefSeq, Jan 2016]





Fibronectin (FN1) Human shRNA Plasmid Kit (Locus ID 2335) - TR304489

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).