

Product datasheet for **TR304772**

EMR2 (ADGRE2) Human shRNA Plasmid Kit (Locus ID 30817)

Product data:

Product Type:	shRNA Plasmids
Product Name:	EMR2 (ADGRE2) Human shRNA Plasmid Kit (Locus ID 30817)
Locus ID:	30817
Synonyms:	CD97; CD312; EMR2; VBU
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ADGRE2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 30817). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001271052 , NM_013447 , NM_152916 , NM_152917 , NM_152918 , NM_152919 , NM_152920 , NM_152921 , NM_013447.1 , NM_013447.2 , NM_013447.3 , NM_152918.1 , NM_001271052.1 , NM_152920.1 , NM_152917.1 , NM_152916.1 , NM_152919.1 , NM_152921.1 , BC127004 , BC034949 , BC053933 , BC082991 , BC127005 , NM_013447.4
UniProt ID:	Q9UHX3
Summary:	This gene encodes a member of the class B seven-span transmembrane (TM7) subfamily of G-protein coupled receptors. These proteins are characterized by an extended extracellular region with a variable number of N-terminal epidermal growth factor-like domains coupled to a TM7 domain via a mucin-like spacer domain. The encoded protein is expressed mainly in myeloid cells where it promotes cell-cell adhesion through interaction with chondroitin sulfate chains. This gene is situated in a cluster of related genes on chromosome 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).