

## Product datasheet for **TR305269**

### CSN7b (COPS7B) Human shRNA Plasmid Kit (Locus ID 64708)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CSN7b (COPS7B) Human shRNA Plasmid Kit (Locus ID 64708)
Locus ID:	64708
Synonyms:	CSN7B; SGN7b
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	COPS7B - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 64708). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001282949</a> , <a href="#">NM_001282950</a> , <a href="#">NM_001282951</a> , <a href="#">NM_001282952</a> , <a href="#">NM_001308381</a> , <a href="#">NM_022730</a> , <a href="#">NM_022730.1</a> , <a href="#">NM_022730.2</a> , <a href="#">NM_022730.3</a> , <a href="#">NM_001282952.1</a> , <a href="#">NM_001282952.2</a> , <a href="#">NM_001282951.1</a> , <a href="#">NM_001282951.2</a> , <a href="#">NM_001282949.1</a> , <a href="#">NM_001282949.2</a> , <a href="#">NM_001282950.1</a> , <a href="#">NM_001282950.2</a> , <a href="#">BC010739</a> , <a href="#">BC091493</a> , <a href="#">NM_001369483</a> , <a href="#">NM_001282949.3</a> , <a href="#">NM_022730.4</a> , <a href="#">NM_001282950.3</a>
UniProt ID:	<a href="#">Q9H9Q2</a>
Summary:	Component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. The complex is also involved in phosphorylation of p53/TP53, JUN, I-kappa-B-alpha/NFKBIA, ITPK1 and IRF8/ICSBP, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).