

## Product datasheet for **TR305372**

### CLCN3 Human shRNA Plasmid Kit (Locus ID 1182)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CLCN3 Human shRNA Plasmid Kit (Locus ID 1182)
Locus ID:	1182
Synonyms:	CIC-3; CLC3
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CLCN3 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 1182). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001243372</a> , <a href="#">NM_001243374</a> , <a href="#">NM_001829</a> , <a href="#">NM_173872</a> , <a href="#">NM_001829.1</a> , <a href="#">NM_001829.2</a> , <a href="#">NM_001829.3</a> , <a href="#">NM_173872.1</a> , <a href="#">NM_173872.2</a> , <a href="#">NM_173872.3</a> , <a href="#">NM_001243372.1</a> , <a href="#">NM_001243374.1</a> , <a href="#">BC026256</a> , <a href="#">BC054877</a> , <a href="#">BC058818</a> , <a href="#">BC136510</a> , <a href="#">BC146697</a> , <a href="#">NM_173872.4</a> , <a href="#">NM_001243372.2</a> , <a href="#">NM_001829.4</a>
UniProt ID:	<a href="#">P51790</a>
Summary:	This gene encodes a member of the voltage-gated chloride channel (CIC) family. The encoded protein is present in all cell types and localized in plasma membranes and in intracellular vesicles. It is a multi-pass membrane protein which contains a CIC domain and two additional C-terminal CBS (cystathionine beta-synthase) domains. The CIC domain catalyzes the selective flow of Cl <sup>-</sup> ions across cell membranes, and the CBS domain may have a regulatory function. This protein plays a role in both acidification and transmitter loading of GABAergic synaptic vesicles, and in smooth muscle cell activation and neointima formation. This protein is required for lysophosphatidic acid (LPA)-activated Cl <sup>-</sup> current activity and fibroblast-to-myofibroblast differentiation. The protein activity is regulated by Ca(2+)/calmodulin-dependent protein kinase II (CaMKII) in glioma cells. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]



[View online »](#)

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).