

Product datasheet for **TR306042**

C1orf135 (AUNIP) Human shRNA Plasmid Kit (Locus ID 79000)

Product data:

Product Type:	shRNA Plasmids
Product Name:	C1orf135 (AUNIP) Human shRNA Plasmid Kit (Locus ID 79000)
Locus ID:	79000
Synonyms:	AIBP; C1orf135
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	AUNIP - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 79000). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC000209 , NM_001287490 , NM_024037 , NM_024037.1 , NM_024037.2 , NM_001287490.1 , BC000209.1 , NM_024037.3
UniProt ID:	Q9H7T9
Summary:	DNA-binding protein that accumulates at DNA double-strand breaks (DSBs) following DNA damage and promotes DNA resection and homologous recombination (PubMed:29042561). Serves as a sensor of DNA damage: binds DNA with a strong preference for DNA substrates that mimic structures generated at stalled replication forks, and anchors RBBP8/CtIP to DSB sites to promote DNA end resection and ensuing homologous recombination repair (PubMed:29042561). Inhibits non-homologous end joining (NHEJ) (PubMed:29042561). Required for the dynamic movement of AURKA at the centrosomes and spindle apparatus during the cell cycle (PubMed:20596670).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).