

## Product datasheet for **TR306428**

### BCCIP Human shRNA Plasmid Kit (Locus ID 56647)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	BCCIP Human shRNA Plasmid Kit (Locus ID 56647)
Locus ID:	56647
Synonyms:	TOK-1; TOK1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	BCCIP - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 56647). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_016567</a> , <a href="#">NM_078468</a> , <a href="#">NM_078469</a> , <a href="#">NM_078468.1</a> , <a href="#">NM_078468.2</a> , <a href="#">NM_078469.1</a> , <a href="#">NM_078469.2</a> , <a href="#">NM_016567.2</a> , <a href="#">NM_016567.3</a> , <a href="#">BC009771</a> , <a href="#">BC009771.1</a> , <a href="#">BM741863</a> , <a href="#">BM905542</a> , <a href="#">BM974496</a> , <a href="#">NM_078469.3</a> , <a href="#">NM_016567.4</a> , <a href="#">NM_078468.3</a>
UniProt ID:	<a href="#">Q9P287</a>
Summary:	This gene product was isolated on the basis of its interaction with BRCA2 and p21 proteins. It is an evolutionarily conserved nuclear protein with multiple interacting domains. The N-terminal half shares moderate homology with regions of calmodulin and M-calpain, suggesting that it may also bind calcium. Functional studies indicate that this protein may be an important cofactor for BRCA2 in tumor suppression, and a modulator of CDK2 kinase activity via p21. This protein has also been implicated in the regulation of BRCA2 and RAD51 nuclear focus formation, double-strand break-induced homologous recombination, and cell cycle progression. Multiple transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).