

## Product datasheet for **TR306527**

### ASXL1 Human shRNA Plasmid Kit (Locus ID 171023)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ASXL1 Human shRNA Plasmid Kit (Locus ID 171023)
Locus ID:	171023
Synonyms:	BOPS; MDS
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ASXL1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 171023). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001164603</a> , <a href="#">NM_015338</a> , <a href="#">NM_015338.2</a> , <a href="#">NM_015338.3</a> , <a href="#">NM_015338.4</a> , <a href="#">NM_015338.5</a> , <a href="#">NM_001164603.1</a> , <a href="#">BC033284</a> , <a href="#">BC064984</a> , <a href="#">BC100280</a> , <a href="#">BC137278</a> , <a href="#">BC137280</a> , <a href="#">NM_001363734</a> , <a href="#">NM_015338.6</a>
UniProt ID:	<a href="#">Q8IXJ9</a>
Summary:	This gene is similar to the Drosophila additional sex combs gene, which encodes a chromatin-binding protein required for normal determination of segment identity in the developing embryo. The protein is a member of the Polycomb group of proteins, which are necessary for the maintenance of stable repression of homeotic and other loci. The protein is thought to disrupt chromatin in localized areas, enhancing transcription of certain genes while repressing the transcription of other genes. The protein encoded by this gene functions as a ligand-dependent co-activator for retinoic acid receptor in cooperation with nuclear receptor coactivator 1. Mutations in this gene are associated with myelodysplastic syndromes and chronic myelomonocytic leukemia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).