

Product datasheet for **TR306622**

ARHGAP11A Human shRNA Plasmid Kit (Locus ID 9824)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ARHGAP11A Human shRNA Plasmid Kit (Locus ID 9824)
Locus ID:	9824
Synonyms:	GAP (1-12)
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ARHGAP11A - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9824). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u>NM_001286479</u> , <u>NM_001286480</u> , <u>NM_014783</u> , <u>NM_199357</u> , <u>NM_014783.1</u> , <u>NM_014783.2</u> , <u>NM_014783.3</u> , <u>NM_014783.4</u> , <u>NM_014783.5</u> , <u>NM_199357.1</u> , <u>NM_199357.2</u> , <u>NM_001286479.1</u> , <u>NM_001286479.2</u> , <u>NM_001286480.1</u> , <u>NM_001286480.2</u> , <u>BC039563</u> , <u>BC039563.2</u> , <u>BC050545</u> , <u>BC063444</u> , <u>NM_001286480.3</u> , <u>NM_014783.6</u>
UniProt ID:	<u>Q6P4F7</u>
Summary:	This gene encodes a member of the Rho GTPase activating protein family. In response to DNA damage, the encoded protein interacts with the p53 tumor suppressor protein and stimulates its tetramerization, which results in cell-cycle arrest and apoptosis. A chromosomal deletion that includes this gene is one cause of Prader-Willi syndrome, and an intronic variant of this gene may be associated with sleep duration in children. This gene is highly expressed in colon cancers and in a human basal-like breast cancer cell line. This gene also produces a ARHGAP11A-SCG5 readthrough transcript and ARHGAP11A-SCG5 protein. [provided by RefSeq, Feb 2019]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).