

## Product datasheet for **TR306658**

### AP1B1 Human shRNA Plasmid Kit (Locus ID 162)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	AP1B1 Human shRNA Plasmid Kit (Locus ID 162)
Locus ID:	162
Synonyms:	ADTB1; AP105A; BAM22; CLAPB2; KIDAR
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	AP1B1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 162). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001127</a> , <a href="#">NM_001166019</a> , <a href="#">NM_145730</a> , <a href="#">NM_001127.1</a> , <a href="#">NM_001127.2</a> , <a href="#">NM_001127.3</a> , <a href="#">NM_145730.1</a> , <a href="#">NM_145730.2</a> , <a href="#">NM_001166019.1</a> , <a href="#">BC046242</a> , <a href="#">BM909530</a> , <a href="#">NM_001127.4</a>
UniProt ID:	<a href="#">Q10567</a>
Summary:	Adaptor protein complex 1 is found at the cytoplasmic face of coated vesicles located at the Golgi complex, where it mediates both the recruitment of clathrin to the membrane and the recognition of sorting signals within the cytosolic tails of transmembrane receptors. This complex is a heterotetramer composed of two large, one medium, and one small adaptin subunit. The protein encoded by this gene serves as one of the large subunits of this complex and is a member of the adaptin protein family. This gene is a candidate meningioma gene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).