

## Product datasheet for **TR306812**

### ADNP Human shRNA Plasmid Kit (Locus ID 23394)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ADNP Human shRNA Plasmid Kit (Locus ID 23394)
Locus ID:	23394
Synonyms:	ADNP1; HVDAS; MRD28
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ADNP - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23394). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001282531</a> , <a href="#">NM_001282532</a> , <a href="#">NM_015339</a> , <a href="#">NM_181442</a> , <a href="#">NM_001347511</a> , <a href="#">NM_015339.1</a> , <a href="#">NM_015339.2</a> , <a href="#">NM_015339.3</a> , <a href="#">NM_015339.4</a> , <a href="#">NM_181442.1</a> , <a href="#">NM_181442.2</a> , <a href="#">NM_181442.3</a> , <a href="#">NM_001282531.1</a> , <a href="#">NM_001282531.2</a> , <a href="#">NM_001282532.1</a> , <a href="#">BC075794</a> , <a href="#">BC075794.1</a> , <a href="#">BC015554</a> , <a href="#">BC090933</a> , <a href="#">BM998017</a> , <a href="#">NM_181442.4</a> , <a href="#">NM_015339.5</a> , <a href="#">NM_001282531.3</a>
UniProt ID:	<a href="#">Q9H2P0</a>
Summary:	Vasoactive intestinal peptide is a neuroprotective factor that has a stimulatory effect on the growth of some tumor cells and an inhibitory effect on others. This gene encodes a protein that is upregulated by vasoactive intestinal peptide and may be involved in its stimulatory effect on certain tumor cells. The encoded protein contains one homeobox and nine zinc finger domains, suggesting that it functions as a transcription factor. This gene is also upregulated in normal proliferative tissues. Finally, the encoded protein may increase the viability of certain cell types through modulation of p53 activity. Alternatively spliced transcript variants encoding the same protein have been described. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).