

## Product datasheet for **TR306833**

### ADAMTS20 Human shRNA Plasmid Kit (Locus ID 80070)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | ADAMTS20 Human shRNA Plasmid Kit (Locus ID 80070)  |
| Locus ID:                 | 80070  |
| Synonyms:                 | ADAM metalloproteinase with thrombospondin type 1 motif, 20; a disintegrin-like and metalloprotease (reprolysin type) with thrombospondin type 1 motif, 20; a disintegrin-like and metalloprotease with thrombospondin type 1 motifs 20; GON-1   |
| Vector:                   | pRS (TR20003)  |
| E. coli Selection:        | Ampicillin   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | ADAMTS20 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 80070). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.  |
| RefSeq:                   | <a href="#">NM_001077205</a> , <a href="#">NM_025003</a> , <a href="#">NM_175851</a> , <a href="#">NM_025003.1</a> , <a href="#">NM_025003.2</a> , <a href="#">NM_025003.3</a> , <a href="#">NM_025003.4</a> , <a href="#">NM_001077205.1</a> , <a href="#">NM_175851.1</a> , <a href="#">NM_025003.5</a>  |
| UniProt ID:               | <a href="#">P59510</a>   |
| Summary:                  | The protein encoded by this gene is a member of the ADAMTS family of zinc-dependent proteases. The encoded protein has a signal peptide that is cleaved to release the mature peptide, which is secreted and found in the extracellular matrix. This protein may be involved in tissue remodeling. [provided by RefSeq, Sep 2011]                  |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> . |


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).