

## Product datasheet for **TR307312**

### **NIP30 (FAM192A) Human shRNA Plasmid Kit (Locus ID 80011)**

#### **Product data:**

Product Type:	shRNA Plasmids
Product Name:	NIP30 (FAM192A) Human shRNA Plasmid Kit (Locus ID 80011)
Locus ID:	80011
Synonyms:	C16orf94; CDA018; CDA10; FAM192A; NIP30; PIP30
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	FAM192A - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 80011). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_024946</a> , <a href="#">NM_001354078</a> , <a href="#">NM_001354079</a> , <a href="#">NM_001354080</a> , <a href="#">NM_001354081</a> , <a href="#">NM_001354082</a> , <a href="#">NM_001354083</a> , <a href="#">NM_001354084</a> , <a href="#">NM_001354085</a> , <a href="#">NM_001354086</a> , <a href="#">NM_001354087</a> , <a href="#">NM_001354088</a> , <a href="#">NM_001354089</a> , <a href="#">NM_001354090</a> , <a href="#">NM_001354091</a> , <a href="#">NM_001354092</a> , <a href="#">NM_001354093</a> , <a href="#">NM_001354094</a> , <a href="#">NM_001354095</a> , <a href="#">NM_001354096</a> , <a href="#">NM_001354097</a> , <a href="#">NM_001354098</a> , <a href="#">NM_001354099</a> , <a href="#">NM_001354100</a> , <a href="#">NM_001354101</a> , <a href="#">NM_001354102</a> , <a href="#">NM_001354103</a> , <a href="#">NM_024946.1</a> , <a href="#">NM_024946.2</a> , <a href="#">BC020536</a> , <a href="#">BC063409</a> , <a href="#">BC071952</a> , <a href="#">NM_024946.4</a>
UniProt ID:	<a href="#">Q9GZU8</a>
Summary:	Promotes the association of the proteasome activator complex subunit PSME3 with the 20S proteasome and regulates its activity. Inhibits PSME3-mediated degradation of some proteasome substrates, probably by affecting their diffusion rate into the catalytic chamber of the proteasome. Also inhibits the interaction of PSME3 with COIL, inhibits accumulation of PSME3 in Cajal bodies and positively regulates the number of Cajal bodies in the nucleus. [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).