

## Product datasheet for **TR307456**

### RAET1L Human shRNA Plasmid Kit (Locus ID 154064)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	RAET1L Human shRNA Plasmid Kit (Locus ID 154064)
Locus ID:	154064
Synonyms:	ULBP6
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RAET1L - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 154064). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_130900</a> , <a href="#">NM_130900.1</a> , <a href="#">NM_130900.2</a> , <a href="#">BC131600</a> , <a href="#">BC148627</a> , <a href="#">BC153191</a> , <a href="#">NM_130900.3</a>
UniProt ID:	<a href="#">Q5VY80</a>
Summary:	RAET1L belongs to the RAET1 family of major histocompatibility complex (MHC) class I-related genes, which are located within a 180-kb cluster on chromosome 6q24.2-q25.3. The REAT1 genes encode glycoproteins that contain extracellular alpha-1 and alpha-2 domains, but they lack the membrane proximal Ig-like alpha-3 domain. Most RAET1 glycoproteins are anchored to the membrane via glycosylphosphatidylinositol (GPI) linkage (Radosavljevic et al., 2002 [PubMed 11827464]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).