

Product datasheet for **TR307683**

TAPT1 Human shRNA Plasmid Kit (Locus ID 202018)

Product data:

Product Type:	shRNA Plasmids
Product Name:	TAPT1 Human shRNA Plasmid Kit (Locus ID 202018)
Locus ID:	202018
Synonyms:	CMVFR; OCLSBG
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TAPT1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 202018). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_153365 , NM_153365.1 , NM_153365.2 , BC066899 , NM_153365.3
UniProt ID:	Q6NXT6
Summary:	This gene encodes a highly conserved protein that localizes to the centrosome and/or ciliary basal body. Mutations in this gene disrupt Golgi morphology and trafficking and normal primary cilium formation and these mutations are congenitally manifested by severe undermineralization of the intra-uterine skeleton. A mutation in the mouse ortholog of this gene results in homeotic, posterior-to-anterior transformations of the axial skeleton which are similar to the phenotype of mouse homeobox C8 gene mutants. In mouse, this gene is thought to function downstream of homeobox C8 to transduce extracellular patterning information during axial skeleton development. [provided by RefSeq, Jan 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).