

## Product datasheet for **TR308376**

### CIAO1 Human shRNA Plasmid Kit (Locus ID 9391)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CIAO1 Human shRNA Plasmid Kit (Locus ID 9391)
Locus ID:	9391
Synonyms:	CIA1; WDR39
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CIAO1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9391). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_004804</a> , <a href="#">NM_004804.1</a> , <a href="#">NM_004804.2</a> , <a href="#">BC032812</a> , <a href="#">BC032812.1</a> , <a href="#">BC001395</a> , <a href="#">NM_004804.3</a>
UniProt ID:	<a href="#">O76071</a>
Summary:	Key component of the cytosolic iron-sulfur protein assembly (CIA) complex, a multiprotein complex that mediates the incorporation of iron-sulfur cluster into extramitochondrial Fe/S proteins (PubMed:17937914, PubMed:23891004). As a CIA complex component, interacts specifically with CIAO2A or CIAO2B and MMS19 to assist different branches of iron-sulfur protein assembly, depending of its interactors. The complex CIAO1:CIAO2B:MMS19 binds to and facilitates the assembly of most cytosolic-nuclear Fe/S proteins. CIAO1:CIAO2A specifically matures ACO1 and stabilizes IREB2 (PubMed:23891004). Seems to specifically modulate the transactivation activity of WT1 (PubMed:9556563). As part of the mitotic spindle-associated MMXD complex it may play a role in chromosome segregation (PubMed:20797633). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).