

## Product datasheet for **TR308401**

### YL1 (VPS72) Human shRNA Plasmid Kit (Locus ID 6944)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	YL1 (VPS72) Human shRNA Plasmid Kit (Locus ID 6944)
Locus ID:	6944
Synonyms:	CFL1; Swc2; TCFL1; YL-1; YL1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	VPS72 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6944). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001271087</a> , <a href="#">NM_001271088</a> , <a href="#">NM_005997</a> , <a href="#">NM_005997.1</a> , <a href="#">NM_005997.2</a> , <a href="#">NM_001271087.1</a> , <a href="#">NM_001271088.1</a> , <a href="#">BC003151</a> , <a href="#">BC003151.1</a> , <a href="#">NM_005997.3</a> , <a href="#">NM_001271088.2</a> , <a href="#">NM_001271087.2</a>
UniProt ID:	<a href="#">Q15906</a>
Summary:	The protein encoded by this gene is a shared subunit of two multi-component complexes, the histone acetyltransferase complex TRRAP/TIP60 as well as the chromatin remodeling SRCAP-containing complex. The TRRAP/TIP60 complex acetylates nucleosomal histones important for transcriptional regulation, double strand DNA break repair and apoptosis. The SRCAP-containing complex catalyzes the exchange of histone H2A with the histone variant Htz1 (H2AFZ) into nucleosomes. This protein may be responsible for binding H2AFZ, which has a role in chromosome segregation. This protein may also have a role in regulating long-term hematopoietic stem cell activity. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Aug 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).