

## Product datasheet for **TR308622**

### TRPM5 Human shRNA Plasmid Kit (Locus ID 29850)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TRPM5 Human shRNA Plasmid Kit (Locus ID 29850)
Locus ID:	29850
Synonyms:	LTRPC5; MTR1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TRPM5 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 29850). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_014555</a> , <a href="#">NM_014555.1</a> , <a href="#">NM_014555.2</a> , <a href="#">NM_014555.3</a> , <a href="#">BC093787</a> , <a href="#">BC093787.1</a> , <a href="#">BC093789</a> , <a href="#">BC143351</a> , <a href="#">BC143353</a>
UniProt ID:	<a href="#">Q9NZQ8</a>
Summary:	This gene encodes a member of the transient receptor potential (TRP) protein family, which is a diverse group of proteins with structural features typical of ion channels. This protein plays an important role in taste transduction, and has characteristics of a calcium-activated, non-selective cation channel that carries Na <sup>+</sup> , K <sup>+</sup> , and Cs <sup>+</sup> ions equally well, but not Ca(2 <sup>+</sup> ) ions. It is activated by lower concentrations of intracellular Ca(2 <sup>+</sup> ), and inhibited by higher concentrations. It is also a highly temperature-sensitive, heat activated channel showing a steep increase of inward currents at temperatures between 15 and 35 degrees Celsius. This gene is located within the Beckwith-Wiedemann syndrome critical region-1 on chromosome 11p15.5, and has been shown to be imprinted, with exclusive expression from the paternal allele. [provided by RefSeq, Oct 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).