

## Product datasheet for **TR30856**

### TGDS Human shRNA Plasmid Kit (Locus ID 23483)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	TGDS Human shRNA Plasmid Kit (Locus ID 23483)
Locus ID:	23483
Synonyms:	CATMANS; SDR2E1; TDPGD
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TGDS - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23483). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001304430</a> , <a href="#">NM_014305</a> , <a href="#">NR_130731</a> , <a href="#">NR_130732</a> , <a href="#">NM_014305.1</a> , <a href="#">NM_014305.2</a> , <a href="#">NM_014305.3</a> , <a href="#">BC033675</a> , <a href="#">BC033675.1</a> , <a href="#">BC005284</a> , <a href="#">NM_014305.4</a>
UniProt ID:	<a href="#">O95455</a>
Summary:	The protein encoded by this gene is a member of the short-chain dehydrogenases/reductases (SDR) superfamily, and is thought to contain a nicotinamide adenine dinucleotide (NAD) binding domain. This large SDR family of enzymes is involved in the metabolism of a variety of compounds, including prostaglandins, retinoids, lipids, steroid hormones, and xenobiotics. Mutations in this gene have been associated with Catel-Manzke syndrome, which is characterized by Pierre Robin sequence, and radial deviation of the index finger due to the presence of an accessory bone between the index finger and its proximal phalanx. Pierre Robin sequence is defined by an undersized jaw, backwards displacement of the tongue base that causes an obstruction of the airways, and can also be associated with a cleft palate. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).