

## Product datasheet for **TR309083**

### SSX2 Human shRNA Plasmid Kit (Locus ID 6757)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SSX2 Human shRNA Plasmid Kit (Locus ID 6757)
Locus ID:	6757
Synonyms:	CT5.2; CT5.2A; HD21; HOM-MEL-40; SSX
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	SSX2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6757). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001278697</a> , <a href="#">NM_003147</a> , <a href="#">NM_175698</a> , <a href="#">NM_003147.1</a> , <a href="#">NM_003147.2</a> , <a href="#">NM_003147.3</a> , <a href="#">NM_003147.4</a> , <a href="#">NM_003147.5</a> , <a href="#">NM_175698.1</a> , <a href="#">NM_175698.2</a> , <a href="#">NM_001278697.1</a> , <a href="#">BC002818</a> , <a href="#">BC002818.2</a> , <a href="#">BC007343</a> , <a href="#">BC007343.2</a> , <a href="#">BC103863</a> , <a href="#">BC016957</a> , <a href="#">BC069313</a> , <a href="#">BC071827</a>
UniProt ID:	<a href="#">Q16385</a>
Summary:	The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy. This gene, and also the SSX1 and SSX4 family members, have been involved in t(X;18)(p11.2;q11.2) translocations that are characteristically found in all synovial sarcomas. This translocation results in the fusion of the synovial sarcoma translocation gene on chromosome 18 to one of the SSX genes on chromosome X. The encoded hybrid proteins are likely responsible for transforming activity. Alternative splicing of this gene results in multiple transcript variants. This gene also has an identical duplicate, GeneID: 727837, located about 45 kb downstream in the opposite orientation on chromosome X. [provided by RefSeq, Jul 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).