

## Product datasheet for **TR309634**

### SCAP2 (SKAP2) Human shRNA Plasmid Kit (Locus ID 8935)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SCAP2 (SKAP2) Human shRNA Plasmid Kit (Locus ID 8935)
Locus ID:	8935
Synonyms:	PRAP; RA70; SAPS; SCAP2; SKAP-HOM; SKAP55R
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	SKAP2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 8935). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001303468</a> , <a href="#">NM_003930</a> , <a href="#">NM_003930.1</a> , <a href="#">NM_003930.2</a> , <a href="#">NM_003930.3</a> , <a href="#">NM_003930.4</a> , <a href="#">BC036044</a> , <a href="#">BC036044.1</a> , <a href="#">BC002893</a> , <a href="#">NM_003930.5</a>
UniProt ID:	<a href="#">O75563</a>
Summary:	The protein encoded by this gene shares homology with Src kinase-associated phosphoprotein 1, and is a substrate of Src family kinases. It is an adaptor protein that is thought to play an essential role in the Src signaling pathway, and in regulating proper activation of the immune system. This protein contains an amino terminal coiled-coil domain for self-dimerization, a pleckstrin homology (PH) domain required for interactions with lipids at the membrane, and a Src homology (SH3) domain at the carboxy terminus. Some reports indicate that this protein inhibits actin polymerization through interactions with actin assembly factors, and might negatively regulate the invasiveness of tumors by modulating actin assembly. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).