

## Product datasheet for **TR309739**

### RPL6 Human shRNA Plasmid Kit (Locus ID 6128)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	RPL6 Human shRNA Plasmid Kit (Locus ID 6128)
Locus ID:	6128
Synonyms:	L6; SHUJUN-2; TAXREB107; TXREB1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RPL6 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6128). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_000970</a> , <a href="#">NM_001024662</a> , <a href="#">NM_001320137</a> , <a href="#">NM_001320138</a> , <a href="#">NM_001320139</a> , <a href="#">NM_001320140</a> , <a href="#">NM_001320141</a> , <a href="#">NM_001320142</a> , <a href="#">NM_001024662.1</a> , <a href="#">NM_001024662.2</a> , <a href="#">NM_000970.1</a> , <a href="#">NM_000970.2</a> , <a href="#">NM_000970.3</a> , <a href="#">NM_000970.4</a> , <a href="#">BC004138</a> , <a href="#">BC004138.2</a> , <a href="#">BC022444</a> , <a href="#">BC022444.1</a> , <a href="#">BC104826</a> , <a href="#">BC104826.1</a> , <a href="#">BC013863</a> , <a href="#">BC020679</a> , <a href="#">BC031009</a> , <a href="#">BC032299</a> , <a href="#">BC070195</a> , <a href="#">BC071912</a> , <a href="#">BC104824</a> , <a href="#">NM_000970.5</a> , <a href="#">NM_001024662.3</a>
UniProt ID:	<a href="#">Q02878</a>
Summary:	This gene encodes a protein component of the 60S ribosomal subunit. This protein can bind specifically to domain C of the tax-responsive enhancer element of human T-cell leukemia virus type 1, and may participate in tax-mediated transactivation of transcription. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).