

## Product datasheet for **TR310452**

### PHKG1 Human shRNA Plasmid Kit (Locus ID 5260)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | PHKG1 Human shRNA Plasmid Kit (Locus ID 5260)   |
| Locus ID:                 | 5260  |
| Synonyms:                 | PHKG  |
| Vector:                   | pRS (TR20003)   |
| E. coli Selection:        | Ampicillin  |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Retroviral plasmids   |
| Components:               | PHKG1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 5260). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.  |
| RefSeq:                   | <a href="#">NM_001258459</a> , <a href="#">NM_001258460</a> , <a href="#">NM_006213</a> , <a href="#">NR_047689</a> , <a href="#">NM_006213.1</a> , <a href="#">NM_006213.2</a> , <a href="#">NM_006213.3</a> , <a href="#">NM_006213.4</a> , <a href="#">NM_001258460.1</a> , <a href="#">NM_001258459.1</a> , <a href="#">BC051327</a> , <a href="#">BC051327.1</a> , <a href="#">BC069655</a> , <a href="#">BC069679</a> , <a href="#">BC069738</a> , <a href="#">BC069754</a> , <a href="#">BC074753</a> , <a href="#">NM_006213.5</a>                |
| UniProt ID:               | <a href="#">Q16816</a>  |
| Summary:                  | This gene is a member of the Ser/Thr protein kinase family and encodes a protein with one protein kinase domain and two calmodulin-binding domains. This protein is the catalytic member of a 16 subunit protein kinase complex which contains equimolar ratios of 4 subunit types. The complex is a crucial glycogenolytic regulatory enzyme. This gene has two pseudogenes at chromosome 7q11.21 and one at chromosome 11p11.12. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).