

Product datasheet for **TR311056**

NUP62 Human shRNA Plasmid Kit (Locus ID 23636)

Product data:

Product Type:	shRNA Plasmids
Product Name:	NUP62 Human shRNA Plasmid Kit (Locus ID 23636)
Locus ID:	23636
Synonyms:	IBSN; p62; SNDI
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	NUP62 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23636). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001193357 , NM_012346 , NM_016553 , NM_153718 , NM_153719 , NM_012346.1 , NM_012346.2 , NM_012346.3 , NM_012346.4 , NM_153718.1 , NM_153718.2 , NM_153718.3 , NM_016553.1 , NM_016553.2 , NM_016553.3 , NM_016553.4 , NM_153719.1 , NM_153719.2 , NM_153719.3 , NM_001193357.1 , BC014842 , BC014842.1 , BC050717 , BC050717.2 , BC095410 , BC095410.1 , BC003663 , BC101104 , BC101105 , BC101106 , BC101107 , BM672281 , NM_001193357.2 , NM_153719.4 , NM_153718.4
UniProt ID:	P37198
Summary:	The nuclear pore complex is a massive structure that extends across the nuclear envelope, forming a gateway that regulates the flow of macromolecules between the nucleus and the cytoplasm. Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells. The protein encoded by this gene is a member of the FG-repeat containing nucleoporins and is localized to the nuclear pore central plug. This protein associates with the importin alpha/beta complex which is involved in the import of proteins containing nuclear localization signals. Multiple transcript variants of this gene encode a single protein isoform. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).