

## Product datasheet for **TR311255**

### CD56 (NCAM1) Human shRNA Plasmid Kit (Locus ID 4684)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CD56 (NCAM1) Human shRNA Plasmid Kit (Locus ID 4684)
Locus ID:	4684
Synonyms:	CD56; MSK39; NCAM
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	NCAM1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 4684). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u><a href="#">NM_000615</a></u> , <u><a href="#">NM_001076682</a></u> , <u><a href="#">NM_001242607</a></u> , <u><a href="#">NM_001242608</a></u> , <u><a href="#">NM_181351</a></u> , <u><a href="#">NM_181351.1</a></u> , <u><a href="#">NM_181351.2</a></u> , <u><a href="#">NM_181351.3</a></u> , <u><a href="#">NM_181351.4</a></u> , <u><a href="#">NM_001076682.1</a></u> , <u><a href="#">NM_001076682.2</a></u> , <u><a href="#">NM_001076682.3</a></u> , <u><a href="#">NM_000615.1</a></u> , <u><a href="#">NM_000615.2</a></u> , <u><a href="#">NM_000615.3</a></u> , <u><a href="#">NM_000615.4</a></u> , <u><a href="#">NM_000615.5</a></u> , <u><a href="#">NM_000615.6</a></u> , <u><a href="#">NM_001242608.1</a></u> , <u><a href="#">NM_001242607.1</a></u> , <u><a href="#">BC047244</a></u> , <u><a href="#">BC047244.1</a></u> , <u><a href="#">BC014205</a></u> , <u><a href="#">BC019845</a></u> , <u><a href="#">BC029119</a></u> , <u><a href="#">NM_001076682.4</a></u> , <u><a href="#">NM_001242608.2</a></u> , <u><a href="#">NM_001242607.2</a></u> , <u><a href="#">NM_000615.7</a></u>
UniProt ID:	<u><a href="#">P13591</a></u>



[View online »](#)

**Summary:**

This gene encodes a cell adhesion protein which is a member of the immunoglobulin superfamily. The encoded protein is involved in cell-to-cell interactions as well as cell-matrix interactions during development and differentiation. The encoded protein plays a role in the development of the nervous system by regulating neurogenesis, neurite outgrowth, and cell migration. This protein is also involved in the expansion of T lymphocytes, B lymphocytes and natural killer (NK) cells which play an important role in immune surveillance. This protein plays a role in signal transduction by interacting with fibroblast growth factor receptors, N-cadherin and other components of the extracellular matrix and by triggering signalling cascades involving FYN-focal adhesion kinase (FAK), mitogen-activated protein kinase (MAPK), and phosphatidylinositol 3-kinase (PI3K). One prominent isoform of this gene, cell surface molecule CD56, plays a role in several myeloproliferative disorders such as acute myeloid leukemia and differential expression of this gene is associated with differential disease progression. For example, increased expression of CD56 is correlated with lower survival in acute myeloid leukemia patients whereas increased severity of COVID-19 is correlated with decreased abundance of CD56-expressing NK cells in peripheral blood. Alternative splicing results in multiple transcript variants encoding distinct protein isoforms. [provided by RefSeq, Aug 2020]

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).