

Product datasheet for **TR311513**

MEF2C Human shRNA Plasmid Kit (Locus ID 4208)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MEF2C Human shRNA Plasmid Kit (Locus ID 4208)
Locus ID:	4208
Synonyms:	C5DELq14.3; DEL5q14.3
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	MEF2C - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 4208). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001131005 , NM_001193347 , NM_001193348 , NM_001193349 , NM_001193350 , NM_001308002 , NM_002397 , NM_002397.1 , NM_002397.2 , NM_002397.3 , NM_002397.4 , NM_001131005.1 , NM_001131005.2 , NM_001193349.1 , NM_001193348.1 , NM_001193350.1 , BC026341 , BC152784 , BC156603 , NM_001363581 , NM_001364332 , NM_001364333 , NM_001364335 , NM_001364337 , NM_001364338 , NM_001364339 , NM_001364340 , NM_001364343 , NM_001364344 , NM_001364346 , NM_001364347 , NM_001364349 , NM_001364353 , NM_001364354 , NM_001364355 , NM_001364329 , NM_001364330 , NM_001364331 , NM_001364334 , NM_001364336 , NM_001364341 , NM_001364342 , NM_001364345 , NM_001364348 , NM_001364350 , NM_001364352 , NM_001364356 , NM_001364357 , NM_001193350.2 , NM_002397.5 , NM_001193349.2
UniProt ID:	Q06413
Summary:	This locus encodes a member of the MADS box transcription enhancer factor 2 (MEF2) family of proteins, which play a role in myogenesis. The encoded protein, MEF2 polypeptide C, has both trans-activating and DNA binding activities. This protein may play a role in maintaining the differentiated state of muscle cells. Mutations and deletions at this locus have been associated with severe cognitive disability, stereotypic movements, epilepsy, and cerebral malformation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2010]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).