

## Product datasheet for **TR311727**

### LILRB5 Human shRNA Plasmid Kit (Locus ID 10990)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	LILRB5 Human shRNA Plasmid Kit (Locus ID 10990)
Locus ID:	10990
Synonyms:	CD85C; LIR-8; LIR8
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	LILRB5 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 10990). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001081442</a> , <a href="#">NM_001081443</a> , <a href="#">NM_001304457</a> , <a href="#">NM_006840</a> , <a href="#">NM_001081442.1</a> , <a href="#">NM_001081442.2</a> , <a href="#">NM_006840.1</a> , <a href="#">NM_006840.2</a> , <a href="#">NM_006840.3</a> , <a href="#">NM_006840.4</a> , <a href="#">NM_001081443.1</a> , <a href="#">NM_001081443.2</a> , <a href="#">BC025704</a> , <a href="#">BC025704.1</a> , <a href="#">NM_001081442.3</a>
UniProt ID:	<a href="#">O75023</a>
Summary:	This gene is a member of the leukocyte immunoglobulin-like receptor (LIR) family, which is found in a gene cluster at chromosomal region 19q13.4. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). Several other LIR subfamily B receptors are expressed on immune cells where they bind to MHC class I molecules on antigen-presenting cells and inhibit stimulation of an immune response. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).