

Product datasheet for **TR312010**

KCNJ15 Human shRNA Plasmid Kit (Locus ID 3772)

Product data:

Product Type:	shRNA Plasmids
Product Name:	KCNJ15 Human shRNA Plasmid Kit (Locus ID 3772)
Locus ID:	3772
Synonyms:	IRKK; KIR1.3; KIR4.2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	KCNJ15 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 3772). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001276435 , NM_001276436 , NM_001276437 , NM_001276438 , NM_001276439 , NM_002243 , NM_170736 , NM_170737 , NM_002243.2 , NM_002243.3 , NM_002243.4 , NM_170737.1 , NM_170737.2 , NM_170736.1 , NM_170736.2 , NM_001276435.1 , NM_001276436.1 , NM_001276437.1 , NM_001276438.1 , NM_001276439.1 , BC013327 , BC013327.1 , NM_001276438.2 , NM_170736.3 , NM_001276437.2 , NM_002243.5 , NM_001276439.2
UniProt ID:	Q99712
Summary:	Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein has a greater tendency to allow potassium to flow into a cell rather than out of a cell. Eight transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).