

Product datasheet for TR312160

IL7 Human shRNA Plasmid Kit (Locus ID 3574)

Product data:

Product Type: shRNA Plasmids

Product Name: IL7 Human shRNA Plasmid Kit (Locus ID 3574)

Locus ID: 3574 Synonyms: IL-7

Vector: pRS (TR20003)

E. coli Selection: Ampicillin

Mammalian Cell Puromycin

Selection:

Format: Retroviral plasmids

Components: IL7 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 3574).

5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.

RefSeq: NM 000880, NM 001199886, NM 001199887, NM 001199888, NM 000880.1, NM 000880.2,

NM 000880.3, NM 001199888.1, NM 001199886.1, NM 001199887.1, BC047698, BC047698.1,

BC032487, NM 000880.4

UniProt ID: P13232

Summary: The protein encoded by this gene is a cytokine important for B and T cell development. This

cytokine and the hepatocyte growth factor (HGF) form a heterodimer that functions as a prepro-B cell growth-stimulating factor. IL7 is found to be a cofactor for V(D)J rearrangement of

the T cell receptor beta (TCRB) during early T cell development. This cytokine can be produced locally by intestinal epithelial and epithelial goblet cells, and may serve as a regulatory factor for intestinal mucosal lymphocytes. IL7 plays an essential role in lymphoid

cell survival, and in the maintenance of naive and memory T cells. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional splice variants have been described but their presence in normal tissues has not been confirmed. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) infection can be a potent inducer of

proinflammatory cytokines and chemokines which may defend against the infection, but may also mediate destructive lung injury. Elevated serum IL7 levels, together with several other circulating cytokines and chemokines, has been found to be associated with the severity of

Coronavirus Disease 19 (COVID-19). [provided by RefSeq, Jul 2020]



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com





shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).